



Med Exchange International, Inc.

MedExchange Inc.

2000 Philips Eclipse 1.5T MRI

2000 Philips Eclipse 1.5T PD250 ISS

Short Bore

Passive shielding

PD250S gradients

Coils: Head/ Atlas, Body Flux ISS, Breast, Small Joint, Large Joint, Shoulder, Knee, Body coil, Flux Body

HP8200 computer VIA5.0 software

Phased Array, Perfusion, Diffusion, Neuro, MRA, FSE, Turbo, EPI

Single Cryogen

Chiller

4GB optical storage devise

Located in CA

[Vendor Information](#)